

ERG Thread

- | | |
|-----------------------------------|-------------------------------|
| •FR SCHEMA DE BRANCHEMENT | •PT ESQUEMA DE CONEXÃO |
| •IT SCHEMA DI COLLEGAMENTO | •EN CONNECTION SCHEME |
| •ES ESQUEMA DE CONEXIÓN | •DE VERBINDUNGSCHEMA |
| | •NL VERBINDUNGSCHEMA |

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ing of electroretinography (ERG), to support the diagnosis of retinal dysfunctions. It can be used to perform either full-field ERG or multifocal ERG or Pattern ERG assessment.

INDICATION FOR USE
 The ERG-Thread supports the diagnosis of retinal dysfunctions, such as the assessment of retinal function impairment that can occur with retinal disease and/or neurological conditions.

PATIENT POPULATION, USER(S) AND USE ENVIRONMENT
 The target patients are adults and adolescents that require an ERG testing. The user(s) profile are ophthalmologists and optometrists (under the supervision of an ophthalmologist) practicing in ophthalmic hospitals, clinics, and medical settings.

CONTRAINDICATIONS
 The ERG-Thread is not indicated for use on patients with cases of trauma or disease of the anterior segment of the eye. In the case of ocular pathology, its use requires the consent of an ophthalmologist.

PRECAUTIONS
DO NOT USE THE DEVICE IF THE PACKAGING IS DAMAGED, IF THE STERILISATION MARKER IS NOT RED OR IF THE DEVICE IS NOT STORED WITHIN THE CORRECT TEMPERATURE RANGE.

USE THE DEVICE ONLY FOR CLINICAL ELECTRORETINOGRAPHY (transient use). Human use of the ERG-Thread requires equipment that complies with all applicable medical device safety standards, to include IEC 60601-1 (please see the connection scheme) *.

Please read the instructions below carefully and follow the safety precautions.

- Due to the risk of irritation of the eye, it should be ensured that are kept in a stable position for the duration of the test. Eyelid retractor may be used if required. Patients are recommended not to blink with the eye during the test.
- The manipulation of the ERG-Thread requires the respect of the rules of the art, in particular the cleanliness in the handling of this sterile object and its unique use.
- Mydriatic drops shall be used to dilate the pupils for better light response of the eye to the stimuli.
- It is suitable to keep the patient in supine position during the test.
- In addition of the usual placement process, a topical anaesthesia is recommended prior testing, but general anesthesia should be avoided.
- This electrode should be placed and removed gently. If pain occurs or persists a few hours after the handling, ophthalmological control is considered necessary.
- Redness and mild irritation of the eye, without exacerbation of the condition, is possible after performing the procedure. In the case of a persisting impairment, an ophthalmologic control is necessary.
- The ERG-Thread is not intended for multiple use. If re-used, an infection of the eye can occur, and the measuring results are no longer guaranteed.

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EN DO NOT RE-USE (infection risk)
 MADE IN CH

DEVICE DESCRIPTION
 The ERG-Thread is an electrode used to transmit the electrical signal elicited by light stimulation of the retina, suitable for human application.

INTENDED USE
 ERG-Thread is an electrode intended to be placed on the corneal or conjunctival fornix or scleral area to transmit electrical signals for subsequent measurement and record-