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MADE IN CH

#### DEVICE DESCRIPTION

The ERG-Thread is an electrode used to transmit the electrical signal elicited by light stimulation of the retina, suitable for human application.

#### INTENDED USE

ERG-Thread is an electrode intended to be placed on the corneal or conjunctival fornix or scleral area to transmit electrical signals for subsequent measurement and record-

ing of electroretinography (ERG), to support the diagnosis of retinal dysfunctions. It can be used to perform either full-field ERG or multifocal ERG or Pattern ERG assessment.

#### INDICATION FOR USE

The ERG-Thread supports the diagnosis of retinal dysfunctions, such as the assessment of retinal function impairment that can occur with retinal disease and/or neurological conditions.

#### PATIENT POPULATION, USER(S) AND USE ENVIRONMENT

The target patients are adults and adolescents that require an ERG testing. The user(s) profile are ophthalmologists and optometrists (under the supervision of an ophthalmologist) practicing in ophthalmic hospitals, clinics, and medical settings.

#### CONTRAINDICATIONS

The ERG-Thread is not indicated for use on patients with cases of trauma or disease of the anterior segment of the eye. In the case of ocular pathology, its use requires the consent of an ophthalmologist.

#### PRECAUTIONS

**DO NOT USE THE DEVICE IF THE PACKAGING IS DAMAGED, IF THE STERILISATION MARKER IS NOT RED OR IF THE DEVICE IS NOT STORED WITHIN THE CORRECT TEMPERATURE RANGE.**

A topical anaesthesia is recommended prior performing the test; Mydriatic drops shall be used to dilate the pupils for better light response of the eye to the stimuli.

1 Open the package.

2 Remove the two sticky patches and affix the end patch containing the connection wire onto the temporal side of the eye. Drape the connection wire over the patient's ear.

Please read the instructions below carefully and follow the safety precautions.

4 Due to the risk of irritation of the eye, it should be ensured that are kept in a stable position for the duration of the test. Eyelid retractor may be used if required. Patients are recommended not to blink with the eye during the test.

5 As the manipulation of the ERG-Thread requires the respect of the rules of the art, in particular the cleanliness in the handling of this sterile object and its unique use.

6 Stick the other sticky pad onto the nose near the eye.

7 Release the eyelid. The fiber should now be either slightly behind the lid or just at the top of the lower lid in contact with the eye. You may need to perform this step again until the fiber sits correctly.

8 Finally plug the connection wire into the montage\*.

9 As ERG-Thread is a single-use device, it must be segregated from other reusable devices. Once used, it should be disposed and not be returned to a decontamination facility for reprocessing.

For any questions, suggestions, or to report any incident concerning ERG-Thread, please get in touch with (and report any incident to the authority having jurisdiction): Fabrinal

- The ERG-Thread is not indicated to be used more than one hour. Feeling of discomfort could be witnessed during the test.

#### SIDE EFFECTS

It may be that a slight discomfort is experienced during the procedure and the product triggers a slight reddening or irritation of the eye. In the case of a persisting impairment, an ophthalmologic control is necessary.

#### EYE DROP TYPES SUITABLE FOR THE PROCEDURE

ERG-Thread shall be used with eye drops administered by a physician. The selection and use of the solution is the responsibility of the practitioner. Before using any solution, users should familiarize themselves with the eye drop manufacturer's instructions and use the solution accordingly.

The compatibility of ERG-Thread has been verified with the following eye drops: 0.4% Oxybuprocaine, and 0.5% Tropicamide. As a general rule, the compatible eye drops are aqueous solutions with a pH in the range 4 to 8. It is recommended to use preservative-free solutions to prevent interactions with the eye.

#### PROCEDURE

The electrode can be inserted into the patient's eye if you follow these instructions:  
 ! A topical anaesthesia is recommended prior performing the test; Mydriatic drops shall be used to dilate the pupils for better light response of the eye to the stimuli.

1 Open the package.

2 Remove the two sticky patches and affix the end patch containing the connection wire onto the temporal side of the eye. Drape the connection wire over the patient's ear.

3 Ask the patient to look up and slightly pull down their lower eyelid.

4 Drape the fiber over the eye just above the lid, making sure as much of the fiber as possible is in contact with the eye.

5 Release the eyelid. The fiber should now be either slightly behind the lid or just at the top of the lower lid in contact with the eye. You may need to perform this step again until the fiber sits correctly.

6 Stick the other sticky pad onto the nose near the eye.

7 Check that the fiber is not too taut. You might ask the patient if the electrode is comfortable (after a few minutes they will hardly notice its presence). If in doubt move the external sticky pad slightly nearer to the eye to reduce tension.

8 Finally plug the connection wire into the montage\*.

9 As ERG-Thread is a single-use device, it must be segregated from other reusable devices. Once used, it should be disposed and not be returned to a decontamination facility for reprocessing.

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