# ERG Jet™


DO NOT RE-USE MADE IN CH

## MODE D'EMPLOI

### MANUALE

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### MANUAL DE USO

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# ERG Jet™

## INSTRUCTION MANUAL BEDIENUNGSANLEITUNG BEDIENINGS- HANDLEIDING

DO NOT RE-USE MADE IN CH  
+25°C CE 0123  
+16°C MD   
STERILE R REF F-06

EN DO NOT RE-USE (infection risk) MADE IN SWITZERLAND

### INTENDED USE

The ERG-Jet is a contact lens electrode intended to be applied on the cornea and to be used for the measurement and recording of electroretinography (ERG) signals, to support the diagnosis of retinal dysfunctions. It can be used to perform either full-field or multifocal ERG assessment.

### INDICATION FOR USE

The ERG-Jet supports the diagnosis of retinal dysfunctions, such as the assessment of retinal function impairment that can occur with retinal disease and/or neurological conditions.

### PATIENT POPULATION, USER(S) AND USE ENVIRONMENT

The target patients are adults and adolescents that require an ERG testing. The user(s) profile are ophthalmologists and optometrists (under the supervision of an ophthalmologist) practicing in ophthalmic hospitals, clinics, and medical settings.

### CONTRAINDICATIONS

ERG-Jet is not intended for use on patients with cases of trauma or disease of the anterior segment of the eye. In the case of ocular and corneal pathologies, its use requires the consent of an ophthalmologist.

### PRECAUTIONS

DO NOT USE THE LENSES IF THE PACKAGING IS DAMAGED, IF THE STERILISATION MARKER IS NOT RED OR IF THE DEVICE IS NOT STORED WITHIN THE CORRECT TEMPERATURE RANGE.

USE THE LENSES ONLY FOR CLINICAL ELECTRORETINOGRAPHY (transient use). Human use of the ERG-Jet requires equipment that complies with all applicable medical device safety standards, to include IEC 60601-1 (please see the connection scheme)\*.

Due to the risk of irritation of the cornea, it should be ensured that patients are kept in a stable position for the duration of the test. Eyelid retractor may be used if required.

Patients are recommended not to blink with the eye during the test.

The manipulation of the ERG-Jet is identical to that of a corneal lens requiring the respect of the rules of the art. In particular, the cleanliness in the handling of this sterile object and the unique use. Mydriatic drops shall be used to dilate the pupils for better light response of the eye to the stimuli.

This lens should be placed and removed gently. If pain occurs or persists a few hours after the handling, ophthalmological control is considered necessary.

Redness and mild irritation of the eye, without exacerbation of the condition, is possible after performing the procedure.

The inner concavity of the contact lens shall be filled with a protective gel prior posing the electrode on the cornea (please refer to the procedure). Blurred vision could be witnessed during the testing when using protective gel.

It is suitable to keep the patient in supine position during the test. In addition of the usual placement process of a conventional contact lens, a topical cornea anesthesia is recommended prior testing, but general anesthesia should be avoided.

The ERG-Jet is not intended for multiple use. If re-used, an infection of the cornea can occur, and the measuring results are no longer guaranteed.

The ERG-Jet is not indicated to be used more than one hour.

Feeling of discomfort could be witnessed during the test.

### EYE DROP TYPES SUITABLE FOR THE PROCEDURE

ERG-Jet shall be used with eye drops administered by a physician.

The selection and use of the solution is the responsibility of the practitioner. Before using any solution, users should familiarize themselves with the eye drop manufacturer's instructions and use the solution accordingly. As a guide, compatibility between the eye drops and commercially available contact lenses can be used based on the eye drop manufacturer's instructions for use.

The compatibility of ERG-Jet has been verified with the following eye drops: 2% Methylcel, 0.4% Oxybuprocaine, and 0.5% Tropicamide. As a general rule, the compatible eye drops are aqueous solutions with a pH in the range 4 to 8. It is recommended to use preservative-free solutions to prevent possible detachment of the lens during measurement and interaction with the eye.

### SIDE EFFECTS

To date, no undesirable side effects have been observed. However, the following side effects could occur:

- Short term blurring of vision during the procedure, without lasting effect

- Short term discomfort during the procedure, without lasting effect

- Short term corneal irritation, without lasting effect

- Short term corneal abrasion, without lasting effect

### PROPERTIES

#### 1-OPTICAL PROPERTIES

Transmission of the complete visible light spectrum


#### 2-MECHANICAL PROPERTIES

Weight = 1.5 g

#### 3-ELECTRICAL PROPERTIES

Resistance = approx. 12 ohms

### PROCEDURE

1 Anesthetize the cornea. Hold the concave surface of the ERG-jet lens upwards and fill it to slightly below the rim with a 2-3% sterile methylcellulose solution.

2 Set ERG-Jet lens onto the cornea by the same procedure used for a conventional contact lens. First, lift the patient's upper lid and have him/her look slightly downwards. Gently rest the lens front edge (with the concave surface upwards) on the lower lid, then tilt the lens backwards onto the cornea.

3 Attach the connection wire to the patient's face using a medical grade adhesive tape.

4 To remove the lens, pull the upper eyelid upwards, then pull the lead wire gently upwards to lift the lower edge of the lens out of the eye.

5 As ERG-Jet is a single-use device, it must be segregated from other reusable devices. Once used, it should be disposed and not be returned to a decontamination facility for reprocessing.

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