INSTRUCTION FOR USE

Paper copy free of charge within 7 calendar days upon request: info@fabrinal.ch

Rev.: G06/26

INTENDED USE

Tran Goniowash® is a single-use irrigation cannula for washing the anterior chamber of the eye, especially the trabecular meshwork and the irido-corneal angle following cataract surgery. Tran Goniowash® is compatible with the Goniowash washing technique.

INDICATION FOR USE

Tran Goniowash® is indicated for use following cataract surgery combined with Goniowash in the eyes of patients with observed pseudoexfoliation syndrome.

CLINICAL BENEFIT TO BE EXPECTED

In patients with pseudoexfoliation syndrome, cataract surgery combined with the use of Tran Goniowash for the Goniowash technique results in up to 40% reduction in intraocular pressure (IOP) from baseline. This clinical performance is maintained for up to 5 years after surgery and translates into a direct clinical benefit for patients, reflected in an 18% to 82% reduction in the use of hypotensive medication.

PATIENT POPULATION, USER(S) AND USE ENVIRONMENT

Patients requiring cataract surgery who suffer from pseudoexfoliation syndrome.

Tran Goniowash® is to be used in a surgical setting by professional

ophthalmologists trained to conduct cataract surgery.

Tran Goniowash® is designed with a female bore fitting male connectors of BSS dispensers irrigation lines.

CONTRAINDICATIONS

The use of Tran Goniowash® is contraindicated for patients with the following conditions:

- Pathology of the uvea: Inflammation, trauma, atrophy, neo- or hyper-vascularization.
- Absolute glaucoma or neovascular glaucoma.
- Ocular herpes affecting the cornea, inflammation of the anterior chamber, of the uvea or of the retina.
- Ocular trauma (e.g.: Contusion, cut, burn, hemorrhage in the anterior chamber).
- Patients who have had an iridectomy or laser iridotomy.

PRECAUTIONS

DO NOT USE THE TRAN GONIOWASH® IF THE PACKAGING IS DAMAGED.

Please read the instructions below carefully and follow the safety precautions.

- In any case, the physician should only operate after having read all the documentation accompanying the product, including instructions for use.

- Only a surgeon specialized in ophthalmology is allowed to use Tran Goniowash®.
- The use of Tran Goniowash® and the results obtained are under the sole responsibility of the physician.
- Only the physician is competent to evaluate the relevant indications and contraindications of the technique to be employed.
- Tran Goniowash® is designed for single use and cannot be reused or resterilised.
- Discard Tran Goniowash® immediately after use in the medical waste container.
- Wearing disposable sterile gloves is strongly recommended to use the device.
- Use the Tran Goniowash® only for washing the anterior chamber of the eye, especially the trabecular meshwork and the irido-corneal angle following cataract surgery (transient use).

Human use of the Tran Goniowash® requires a phacoemulsification device that complies with all applicable medical device safety standards, to include IEC 60601-1.

PROCEDURE

- 1** At the end of cataract surgery, connect the Tran Goniowash® at the irrigation part of the phacoemulsifier. The BSS should be placed at an average height of 80 cm from the eyeball. Introduce the tip of the cannula into the anterior chamber at the corneoscleral junction.
- 2** Perform the washing of the anterior chamber of the eye being careful not to come in direct contact with the tissues of the inside of the eye.
- 3** As Tran Goniowash is a single-use device, it must be segregated from reusable devices. Once used, it should be disposed and not be returned to a

decontamination facility for reprocessing.

SIDE EFFECTS

In case of contact between the Tran Goniowash® and tissues of the interior of the eye the following complications may occur:

- Temporary deformation of the pupil.
- Thinning, atrophy, inflammation or perforation of the iris.
- Iridodialysis.
- Minimal and brief haemorrhage in the anterior chamber (blood reflux and blood in the periphery of the iris).
- Inflammation of the cornea.
- Depigmentation, pigment dispersion.
- Tearing and detachment of the corneal endothelium.

Moreover, there are also the major complications following cataract surgery:

- Retinal detachment.
- IOL dislocation.
- Corneal decompensation.
- Diplopia.
- Blindness.
- Infections (1 to 3 cases out of 1000).
- Eye trauma caused by the patient or their relations.
- Scar insufficiently waterproof.
- Partial ptosis of the upper eyelid.
- Subconjunctival or eyelid haematoma.
- Perceiving flying spots.
- Ocular inflammation.
- Increase in intraocular pressure.
- Astigmatism.


Other complications following cataract surgery:

- Iris trauma.
- Descemet membrane tears or detachment.

- Posterior capsular tear or zonular rupture.

- Retained lens fragments.

For any questions, suggestions, or to report any incident concerning Tran Goniowash®, please get in touch with (and report any incident to the authority having jurisdiction): Fabrinal

VERANEX GERMANY GMBH
LANDSBERGER STR. 302, 80687 MUNICH
GERMANY SRN : DE-AR-000005578

FABRINAL SA, Rue de la Tuilerie 42
CH-2300 La Chaux-de-Fonds, Switzerland
+41 32 968 09 08, info@fabrinal.ch

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